

Anti-Human/Mouse/Rat RPH3A Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05680

3th Edition

Description:Protein transport. Probably involved with Ras-related protein Rab-3A in synaptic vesicle traffic and/or synaptic vesicle fusion. Could play a role in neurotransmitter release by regulating membrane flow in the nerve terminal.

Antigen:Synthesized peptide derived from human Rabphilin-3A around the non-phosphorylation site of Ser237.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: