

Anti-Human RFTN1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05701

3th Edition

Description: May play a pivotal role in the formation and/or maintenance of lipid rafts. May regulate B-cell antigen receptor-mediated signaling.

Antigen: Synthesized peptide derived from the C-terminal region of human Raftlin

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Specifically expressed by B-cells (at protein level).