

Anti-Human/Mouse/Rat FRK Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05706

3th Edition

Description: Positively regulates PTEN protein stability through phosphorylation of PTEN on 'Tyr-336', which in turn prevents its ubiquitination and degradation possibly by reducing its binding to NEDD4. May function as a tumor suppressor.

Antigen: Synthesized peptide derived from the Internal region of human Rak

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Restricted to cells lines derived from tissues of lymphoid, brain, breast, colon and bladder origin.