

Anti-Human/Rat RASSF4 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05714

3th Edition

Description: The function of this gene has not yet been determined but may involve a role in tumor suppression. Alternative splicing of this gene results in several transcript variants; however, most of the variants have not been fully described. RASSF4 (Ras Association Domain Family Member 4) is a Protein Coding gene. Among its related pathways are Hippo signaling pathway -multiple species. An important paralog of this gene is RASSF2.

Antigen: Synthesized peptide derived from the Internal region of human RASSF4

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: