



**Anti-Human/Mouse/Rat RBM16 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05722**

---

3th Edition

**Description:** SCAF8 (SR-Related CTD Associated Factor 8) is a Protein Coding gene. GO annotations related to this gene include nucleic acid binding and nucleotide binding. An important paralog of this gene is SCAF4.

**Antigen:** Synthesized peptide derived from the Internal region of human RBM16

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**