

**Anti-Human/Mouse/Rat RBM34 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05723**

---

3th Edition

**Description:** This gene encodes a member of the RNA-binding motif family of RNA recognition motif proteins. The encoded protein contains an RNA-binding domain made up of two RNA recognition motif subdomains referred to as RNA recognition motif-1 and RNA recognition motif-2. Alternative splicing results in multiple transcript variants.

**Antigen:** Synthesized peptide derived from the C-terminal region of human RBM34

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**