

**Anti-Human/Mouse/Rat TMBIM1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05728**

---

3th Edition

**Description:**TMBIM1 (Transmembrane BAX Inhibitor Motif Containing 1) is a Protein Coding gene. Among its related pathways are Innate Immune System. GO annotations related to this gene include death receptor binding. An important paralog of this gene is FAIM2.

**Antigen:**Synthesized peptide derived from the C-terminal region of human RECS1

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**