

Anti-Human RPTN Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05732

3th Edition

Description:RPTN (Repetin) is a Protein Coding gene. Diseases associated with RPTN include Autosomal Recessive Congenital Ichthyosis. Among its related pathways are Keratinization and Developmental Biology. GO annotations related to this gene include calcium ion binding. An important paralog of this gene is FLG2.

Antigen:Synthesized peptide derived from the C-terminal region of human Repetin

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expression is scattered in the normal epidermis but strong in the acrosyngium, the inner hair root sheath and in the filiform papilli of the tongue.