

**Anti-Human/Mouse RHOBTB3 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05748**

---

3th Edition

**Description:** RHOBTB3 is a member of the evolutionarily conserved RHOBTB subfamily of Rho GTPases. For background information on RHOBTBs, see RHOBTB1 (MIM 607351). RHOBTB3 (Rho Related BTB Domain Containing 3) is a Protein Coding gene. Among its related pathways are Vesicle-mediated transport and Retrograde transport at the Trans-Golgi-Network. GO annotations related to this gene include GTP binding and Rab GTPase binding.

**Antigen:** Synthesized peptide derived from the Internal region of human RHOBTB3

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Ubiquitous. Highly expressed in neural and cardiac tissues, pancreas, placenta and testis.