

Anti-Human RHPN1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05751

3th Edition

Description:RHPN1 (Rhopilin Rho GTPase Binding Protein 1) is a Protein Coding gene. Among its related pathways are Ectoderm Differentiation and Actin Nucleation by ARP-WASP Complex. An important paralog of this gene is RHPN2.

Antigen:Synthesized peptide derived from the Internal region of human Rhophilin

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: