

Anti-Human RPAP2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05790

3th Edition

Description:RPAP2 (RNA Polymerase II Associated Protein 2) is a Protein Coding gene. Among its related pathways are Gene Expression and RNA polymerase II transcribes snRNA genes. GO annotations related to this gene include CTD phosphatase activity.

Antigen:Synthesized peptide derived from the C-terminal region of human RPAP2

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: