

Anti-Human/Mouse RASL10A Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05792

3th Edition

Description:RASL10A (RAS Like Family 10 Member A) is a Protein Coding gene. GO annotations related to this gene include GTP binding and GTPase activity. An important paralog of this gene is RASL10B.

Antigen:Synthesized peptide derived from the Internal region of human RRP22

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expression appears to be strictly limited to the central nervous system.