

**Anti-Human/Mouse RSAD1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05794**

---

3th Edition

**Description:**RSAD1 (Radical S-Adenosyl Methionine Domain Containing 1) is a Protein Coding gene. GO annotations related to this gene include 4 iron, 4 sulfur cluster binding and coproporphyrinogen oxidase activity.

**Antigen:**Synthesized peptide derived from the Internal region of human RSAD1

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**