

**Anti-Human S100A7L2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05800**

---

3th Edition

**Description:**S100A7L2 (S100 Calcium Binding Protein A7 Like 2) is a Protein Coding gene. GO annotations related to this gene include calcium ion binding. An important paralog of this gene is S100A7A.

**Antigen:**Synthesized peptide derived from the Internal region of human S-100A7L2

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**