

**Anti-Human S100Z Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05801**

---

3th Edition

**Description:**Members of the S100 protein family contain 2 calcium-binding EF-hands and exhibit cell-type specific expression patterns. For additional background information on S100 proteins, see MIM 114085.S100Z (S100 Calcium Binding Protein Z) is a Protein Coding gene. GO annotations related to this gene include calcium ion binding. An important paralog of this gene is S100A1.

**Antigen:**Synthesized peptide derived from the N-terminal region of human S-100Z

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**