

Anti-Human/Mouse/Rat DLGAP1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05807

3th Edition

Description:DLGAP1 (DLG Associated Protein 1) is a Protein Coding gene. Diseases associated with DLGAP1 include Retinitis Pigmentosa 55 and Obsessive-Compulsive Disorder. Among its related pathways are Protein-protein interactions at synapses and Neuroscience. An important paralog of this gene is DLGAP2.

Antigen:Synthesized peptide derived from the C-terminal region of human SAPAP1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed in brain.