

Anti-Human/Mouse SCAF1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05809

3th Edition

Description:SCAF1 (SR-Related CTD Associated Factor 1) is a Protein Coding gene. GO annotations related to this gene include RNA binding and protein domain specific binding. An important paralog of this gene is PHRF1.

Antigen:Synthesized peptide derived from the N-terminal region of human SCAF1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Ubiquitous. Highly expressed in fetal brain and liver, poorly expressed in salivary gland, heart, skin and ovary. Expressed in colorectal carcinomas and ovarian cancers. Overexpressed in colorectal carcinomas as compared to normal colonic mucosa.