

**Anti-Human/Mouse/Rat SCFD1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05812**

---

3th Edition

**Description:**SCFD1 (Sec1 Family Domain Containing 1) is a Protein Coding gene. Among its related pathways are Vesicle-mediated transport and Transport to the Golgi and subsequent modification. GO annotations related to this gene include protein N-terminus binding and syntaxin binding.

**Antigen:**Synthesized peptide derived from the C-terminal region of human SCFD1

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**