

**Anti-Human/Mouse/Rat SDHA Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05815**

---

3th Edition

**Description:** Flavoprotein (FP) subunit of succinate dehydrogenase (SDH) that is involved in complex II of the mitochondrial electron transport chain and is responsible for transferring electrons from succinate to ubiquinone (coenzyme Q) (PubMed:24781757). Can act as a

**Antigen:** Synthesized peptide derived from the C-terminal region of human SDHA.

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**