

Anti-Human/Mouse/Rat SENP8 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05824

3th Edition

Description: This gene encodes a cysteine protease that is a member of the sentrin-specific protease family. The encoded protein is involved in processing and deconjugation of the ubiquitin-like protein termed, neural precursor cell expressed developmentally downregulated 8. Alternate splicing results in multiple transcript variants.

Antigen: Synthesized peptide derived from the C-terminal region of human SENP8

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Broadly expressed, with highest levels in kidney and pancreas.