

**Anti-Human SERPINB12 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05834**

---

3th Edition

**Description:**SERPINB12 (Serp Family B Member 12) is a Protein Coding gene. Among its related pathways are Innate Immune System. GO annotations related to this gene include enzyme binding and serine-type endopeptidase inhibitor activity. An important paralog of this gene is SERPINB4.

**Antigen:**Synthesized peptide derived from the Internal region of human Serpin B12

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expressed in many tissues, including brain, bone marrow, lymph node, heart, lung, liver, pancreas, testis, ovary, and intestine.