

Anti-Human/Mouse/Rat SIAH1/2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05850

3th Edition

Description:SIAH1 (Siah E3 Ubiquitin Protein Ligase 1) is a Protein Coding gene. Diseases associated with SIAH1 include Parkinson Disease, Late-Onset. Among its related pathways are Proteolysis Role of Parkin in the Ubiquitin-Proteasomal Pathway and Innate Immune System. GO annotations related to this gene include identical protein binding and ubiquitin-protein transferase activity. An important paralog of this gene is SIAH2. SIAH2 (Siah E3 Ubiquitin Protein Ligase 2) is a Protein Coding gene. Among its related pathways are Proteolysis Role of Parkin in the Ubiquitin-Proteasomal Pathway and Innate Immune System. GO annotations related to this gene include ligase activity and transcription corepressor activity. An important paralog of this gene is SIAH1.

Antigen:Synthesized peptide derived from the Internal region of human Siah-1/2

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:SIAH1 is Widely expressed at a low level. Down-regulated in advanced hepatocellular carcinomas.SIAH2 isWidely expressed at low level.