

Anti-Human/Mouse SLC17A2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05860

3th Edition

Description: SLC17A2 (Solute Carrier Family 17 Member 2) is a Protein Coding gene. Among its related pathways are Osteoblast Signaling. GO annotations related to this gene include sodium:phosphate symporter activity. An important paralog of this gene is SLC17A4.

Antigen: Synthesized peptide derived from the Internal region of human SLC17A2

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in the small intestine, kidney, spleen and testis. Not detected in fetal brain, bone marrow, and mammary gland.