

Anti-Human/Mouse/Rat SLC5A3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05879

3th Edition

Description:Prevents intracellular accumulation of high concentrations of myo-inositol (an osmolyte) that result in impairment of cellular function.

Antigen:Synthesized peptide derived from the Internal region of human SMIT

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: