

Anti-Human SUPT3H Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05899

3th Edition

Description:SUPT3H (SPT3 Homolog, SAGA And STAGA Complex Component) is a Protein Coding gene. Diseases associated with SUPT3H include Blastic Plasmacytoid Dendritic Cell. Among its related pathways are Validated targets of C-MYC transcriptional activation and C-MYC pathway. GO annotations related to this gene include protein heterodimerization activity and histone acetyltransferase activity.

Antigen:Synthesized peptide derived from the Internal region of human SPT3

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed in all tissues tested including pancreas, kidney, skeletal muscle, liver, lung, placenta, brain and heart.