

**Anti-Human HTR1E Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05901**

---

3th Edition

**Description:**HTR1E (5-Hydroxytryptamine Receptor 1E) is a Protein Coding gene. Diseases associated with HTR1E include Attention Deficit-Hyperactivity Disorder. Among its related pathways are Serotonergic synapse and Peptide ligand-binding receptors. GO annotations related to this gene include G-protein coupled receptor activity and serotonin binding. An important paralog of this gene is HTR1F.

**Antigen:**Synthesized peptide derived from the Internal region of human SR-1E

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Detected in brain.