

Anti-Human/Mouse/Rat/Monkey HTR1F Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05902

3th Edition

Description:HTR1F (5-Hydroxytryptamine Receptor 1F) is a Protein Coding gene. Diseases associated with HTR1F include Migraine With Or Without Aura 1. Among its related pathways are Serotonergic synapse and Peptide ligand-binding receptors. GO annotations related to this gene include G-protein coupled receptor activity and serotonin binding. An important paralog of this gene is HTR1E.

Antigen:Synthesized peptide derived from the Internal region of human SR-1F

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: