



Anti-Human/Mouse SRSF4 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05916

3th Edition

Description: This gene encodes a member of the arginine/serine-rich splicing factor family. The encoded protein likely functions in mRNA processing.

Antigen: Synthesized peptide derived from the Internal region of human SRp75

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: