

**Anti-Human/Mouse/Rat STAT5A Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05934**

---

3th Edition

**Description:** Carries out a dual function: signal transduction and activation of transcription. Binds to the GAS element and activates PRL-induced transcription.

**Antigen:** Synthesized peptide derived from human Stat5a around the non-phosphorylation site of S780.

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**