

Anti-Human ELP2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05939

3th Edition

Description:Regulates the ligand-dependent activation of STAT3.Acts as subunit of the RNA polymerase II elongator complex, which is a histone acetyltransferase component of the RNA polymerase II (Pol II) holoenzyme and is involved in transcriptional elongation. Elongator may play a role in chromatin remodeling and is involved in acetylation of histones H3 and probably H4.

Antigen:Synthesized peptide derived from the C-terminal region of human Stlp1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: