



Anti-Human/Mouse/Rat MMP3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05943

3th Edition

Description: Can degrade fibronectin, laminin, gelatins of type I, III, IV, and V; collagens III, IV, X, and IX, and cartilage proteoglycans. Activates procollagenase.

Antigen: Synthesized peptide derived from the C-terminal region of human MMP-3

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: