

Anti-Human/Mouse/Rat SYT10 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05954

3th Edition

Description: May be involved in Ca(2+)-dependent exocytosis of secretory vesicles through Ca(2+) and phospholipid binding to the C2 domain or may serve as Ca(2+) sensors in the process of vesicular trafficking and exocytosis.

Antigen: Synthesized peptide derived from the C-terminal region of human Synaptotagmin X

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed only in pancreas, lung and kidney.