

**Anti-Human/Mouse/Rat TRAF3IP3 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05964**

---

3th Edition

**Description:** The gene encodes a protein that mediates cell growth by modulating the c-Jun N-terminal kinase signal transduction pathway. The encoded protein may also interact with a large multiprotein assembly containing the phosphatase 2A catalytic subunit.

**Antigen:** Synthesized peptide derived from the Internal region of human T3JAM

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**