

**Anti-Human TAAR3P Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05965**

---

3th Edition

**Description:**TAAR3P (Trace Amine Associated Receptor 3, Pseudogene) is a Pseudogene. Among its related pathways are Peptide ligand-binding receptors and Signaling by GPCR.

**Antigen:**Synthesized peptide derived from the Internal region of human TAAR3

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Not expressed in the pons, thalamus, globus pallidus, caudate, putamen or cerebellum.