

**Anti-Human/Mouse TNKS1BP1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA05966**

---

3th Edition

**Description:** TNKS1BP1 (Tankyrase 1 Binding Protein 1) is a Protein Coding gene. Among its related pathways are Gene Expression and Deadenylation-dependent mRNA decay. GO annotations related to this gene include enzyme binding and ankyrin binding.

**Antigen:** Synthesized peptide derived from the C-terminal region of human TAB182

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Detected in testis, ovary, lung, skeletal muscle, heart, prostate and pancreas, and at very low levels in brain and peripheral blood leukocytes.