

Anti-Human TBC1D3A/B/C Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05981

3th Edition

Description: Acts as a GTPase activating protein for RAB5. Does not act on RAB4 or RAB11.

Antigen: Synthesized peptide derived from the C-terminal region of human TBC1D3A/B/C

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in liver, skeletal muscle, kidney, pancreas, spleen, testis, ovary, small intestine and peripheral blood leukocytes. Overexpressed in prostate cancers.