

Anti-Human/Mouse TBX15/18 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA05985

3th Edition

Description: TBX15 (T-Box 15) is a Protein Coding gene. Diseases associated with TBX15 include Cousin Syndrome and Acromegaloid Facial Appearance Syndrome. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and protein heterodimerization activity. An important paralog of this gene is TBX18. TBX18 (T-Box 18) is a Protein Coding gene. Diseases associated with TBX18 include Pelviureteric Junction Obstruction and Huntington Disease. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and protein heterodimerization activity. An important paralog of this gene is TBX15.

Antigen: Synthesized peptide derived from the Internal region of human TBX15/18

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: