

Anti-Human/Mouse TNS1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06004

3th Edition

Description:The protein encoded by this gene localizes to focal adhesions, regions of the plasma membrane where the cell attaches to the extracellular matrix. This protein crosslinks actin filaments and contains a Src homology 2 (SH2) domain, which is often found in molecules involved in signal transduction. This protein is a substrate of calpain II. Alternative splicing results in multiple transcript variants encoding different isoforms. TNS1 (Tensin 1) is a Protein Coding gene. Diseases associated with TNS1 include Bannayan-Riley-Ruvalcaba Syndrome and Ruvalcaba Syndrome. Among its related pathways are Focal Adhesion and Integrin Pathway. GO annotations related to this gene include poly(A) RNA binding and actin binding. An important paralog of this gene is TNS3.

Antigen:Synthesized peptide derived from human Tensin-1 around the non-phosphorylation site of Tyr1326.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Ubiquitous.