

Anti-Human/Mouse/Rat TESK2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06007

3th Edition

Description:Dual specificity protein kinase activity catalyzing autophosphorylation and phosphorylation of exogenous substrates on both serine/threonine and tyrosine residues. Phosphorylates cofilin at 'Ser-3'. May play an important role in spermatogenesis.

Antigen:Synthesized peptide derived from the Internal region of human TESK2

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Predominantly expressed in testis and prostate. Found predominantly in non-germinal Sertoli cells.