

Anti-Human TIGD3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06027

3th Edition

Description:The protein encoded by this gene belongs to the tigger subfamily of the pogo superfamily of DNA-mediated transposons in humans. These proteins are related to DNA transposons found in fungi and nematodes, and more distantly to the Tc1 and mariner transposases. They are also very similar to the major mammalian centromere protein B. The exact function of this gene is not known. TIGD3 (Tigger Transposable Element Derived 3) is a Protein Coding gene. GO annotations related to this gene include nucleic acid binding. An important paralog of this gene is TIGD4.

Antigen:Synthesized peptide derived from the C-terminal region of human TIGD3

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: