

Anti-Human/Mouse/Rat TLE4 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06034

3th Edition

Description: TLE4 (Transducin-like enhancer protein 4) is a member of the Groucho family of transcriptional corepressors and inhibits the transcriptional activation mediated by Pax5 and by beta-catenin and TCF family members in Wnt signaling.

Antigen: Synthesized peptide derived from the Internal region of human TLE4

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: In all tissues examined, mostly in brain, and muscle.