

**Anti-Human/Mouse/Rat TNFAIP1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA06040**

---

3th Edition

**Description:** This gene was identified as a gene whose expression can be induced by the tumor necrosis factor alpha (TNF) in umbilical vein endothelial cells. Studies of a similar gene in mouse suggest that the expression of this gene is developmentally regulated in a tissue-specific manner.

**Antigen:** Synthesized peptide derived from the Internal region of human TNF-IP 1

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**