



Anti-Human TP53INP1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06045

3th Edition

Description:In response to double-strand DNA breaks, promotes p53/TP53 phosphorylation on 'Ser-46' and subsequent apoptosis.

Antigen:Synthesized peptide derived from the C-terminal region of human TP53INP1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Ubiquitously expressed.