

**Anti-Human/Mouse/Rat TACSTD2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA06066**

---

3th Edition

**Description:** This intronless gene encodes a carcinoma-associated antigen. This antigen is a cell surface receptor that transduces calcium signals. Mutations of this gene have been associated with gelatinous drop-like corneal dystrophy.

**Antigen:** Synthesized peptide derived from the Internal region of human TROP-2

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Placenta, pancreatic carcinoma cell lines.