

Anti-Mouse/Rat TNNI3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06067

3th Edition

Description: Troponin I is the inhibitory subunit of troponin, the thin filament regulatory complex which confers calcium-sensitivity to striated muscle actomyosin ATPase activity.

Antigen: Synthesized peptide derived from human Troponin I-C around the non-phosphorylation site of S43.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: