

**Anti-Human/Mouse DCT Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA06069**

---

3th Edition

**Description:** Dopachrome tautomerase (dopachrome delta-isomerase, tyrosine-related protein 2), also known as DCT, is a human gene. The human TYRP2 protein has 83% identity and 90% similarity to the mouse sequence and has all the structural characteristics of the tyrosinase protein family, including a signal peptide, 15 conserved cysteine residues, 2 copper-binding domains, and a C-terminal membrane-spanning region. Northern blot analysis showed that TYRP2 is expressed at high levels in human melanoma cells. Its expression is regulated by the microphthalmia-associated transcription factor (MITF).

**Antigen:** Synthesized peptide derived from the Internal region of human TRP2

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**