

**Anti-Human/Mouse TUSC2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA06075**

---

3th Edition

**Description:** This gene is a highly conserved lung cancer candidate gene. No other information about this gene is currently available.

**Antigen:** Synthesized peptide derived from the Internal region of human TUSC2

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Strong expression in heart, lung, skeletal muscle, kidney, and pancreas, followed by brain and liver, lowest levels in placenta.