

Anti-Human/Mouse UBE2T Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06086

3th Edition

Description:The protein encoded by this gene catalyzes the covalent attachment of ubiquitin to protein substrates. Defects in this gene have been associated with Fanconi anemia of complementation group T. Two transcript variants encoding different isoforms have been found for this gene.

Antigen:Synthesized peptide derived from the Internal region of human UBE2T

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: