

Anti-Human UBE3C Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06087

3th Edition

Description:UBE3C (Ubiquitin Protein Ligase E3C) is a Protein Coding gene. Among its related pathways are Innate Immune System and Class I MHC mediated antigen processing and presentation. GO annotations related to this gene include ligase activity and ubiquitin-protein transferase activity. An important paralog of this gene is UBE3B.

Antigen:Synthesized peptide derived from the Internal region of human UBE3C.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Highly expressed in skeletal muscle. Detected at much lower levels in kidney and pancreas.