

**Anti-Human VN1R2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA06108**

---

3th Edition

**Description:** VN1R2 (Vomeronasal 1 Receptor 2) is a Protein Coding gene. GO annotations related to this gene include G-protein coupled receptor activity and pheromone receptor activity. An important paralog of this gene is VN1R4.

**Antigen:** Synthesized peptide derived from the Internal region of human V1RL2

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**